





Recycling Passport

Product name:	OPTIMUS CXA 65/80 1T 400V		
Identification code	9890-000-70381; 9890-000-70385		
Order No.			
Device Type	B2B (Business to business)		
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	http://www.healthcare.philips.com/main/about/sustainability/recycling/xraygenerators.wpd	
Total weight (in Kg)			216
thereof			

steel, iron	iron, low alloy (<5%)	108,752	KG
	iron, high alloy (>5%)	0,142	KG
steel, iron		108,894	KG
nonferrous metals and alloys	copper, -alloy	3,959	KG
	aluminum, -alloy	0,752	KG
	zinc, -alloy	0,438	KG
nonferrous metals and alloys		5,149	KG
glass / ceramics	glass	0,288	KG
	glass / ceramics		0,288 KG
plastics / organic substances	oil	50,069	KG
	thermoplastic	4,976	KG
	thermoset	1,861	KG
	elastomer	0,221	KG
plastics / organic substances		57,127	KG
standard parts	other electronic powered devic	24,988	KG
	printed circuit boards	10,121	KG
	cables	6,614	KG
	mounting parts, attaching part	2,541	KG
	fans	0,432	KG
	motors, pumps	0,276	KG
standard parts		44,972	KG
		216,43	KG

Additional weight for Systemcomponents see page three.

Recycling Passport

Recycle info		
General Hazards (Desinfection)	<p>For dismantling activities Treatment Facilities must consider the national requirements.</p> <p>For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account</p>	
Hazardous Substances  To be removed	<p style="text-align: center;">Substances:</p>	<p style="text-align: center;">Location</p>
Batteries  To be removed	<p style="text-align: center;">Type:</p> <p>Lithium chrome cell, 3V</p> <p>Lithium chrome cell, 3V</p>	<p style="text-align: center;">Location</p> <p>Figure 1 (1)</p> <p>Figure 1 (2)</p>
Special attention 	<p style="text-align: center;">Item</p> <p>Pb is present in the soldering process of PCBs</p>	<p style="text-align: center;">Location</p> <p>Figure 1, 2</p>
Fluids / Gases 	<p style="text-align: center;">Item</p> <p>Transformer Oil, Shell Diala</p> <p>HV-transformer oil contains no PCBs</p>	<p style="text-align: center;">Location</p> <p>Figure 1 (3)</p>

Recycling Passport




Figure 1

- 4.** Room in/out Code Nr. 4512-201-0417. overall weight 1,9 kg. Materials see Recycling Passport.
4.1 Cableset cabinet Code Nr. 4512-201-0567. overall weight 0,65 kg. Materials see Recycling passport.
4.2 System Power distribution Unit Code Nr. 4512-201-0476. overall weight 4,8 kg. Materials see Recycling Passport.