





Recycling Passport

Product name:	OPTIMUS CXA 65/80 1T 480V		
Identification code	9890-000-70391; 9890-000-70395		
Order No.			
Device Type	B2B (Business to business)		
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	http://www.medical.philips.com/main/company/sustainability/recycling/xray_generators/	
Total weight (in Kg) thereof		254	

steel, iron	iron, low alloy (<5%)	110,18	KG
	iron, high alloy (>5%)	0,162	KG
steel, iron		110,342	KG
nonferrous metals and alloys	copper, -alloy	3,959	KG
	aluminum, -alloy	0,752	KG
	zinc, -alloy	0,438	KG
nonferrous metals and alloys		5,149	KG
glass / ceramics	glass	0,288	KG
glass / ceramics		0,288	KG
plastics / organic substances	oil	50,069	KG
	thermoplastic	4,987	KG
	thermoset	3,947	KG
	elastomer	0,221	KG
plastics / organic substances		59,224	KG
standard parts	other electronic powered devic	58,698	KG
	printed circuit boards	10,121	KG
	cables	6,946	KG
	mounting parts, attaching part	2,683	KG
	motors, pumps	0,276	KG
	fans	0,064	KG
standard parts		78,788	KG
		253,791	KG

Additional weight for Systemcomponents see page three.

Recycling Passport

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances  To be removed	Substances:	Location
Batteries  To be removed	Type:	Location
	Lithium chrome cell, 3V	Figure 1 (1)
	Lithium chrome cell, 3V	Figure 1 (2)
Special attention 	Item	Location
	Pb is present in the soldering process of PCBs	Figure 1, 2
Fluids / Gases 	Item	Location
	Transformer Oil, Shell Diala	Figure 1 (3)
	HV-transformer oil contains no PCBs	

Recycling Passport




Figure 1

4. Room in/out Code Nr. 4512-201-0417. overall weight 1,9 kg. Materials see Recycling Passport.

4.1 Cableset cabinet Code Nr. 4512-201-0567. overall weight 0,65 kg. Materials see Recycling passport.

4.2 System Power distribution Unit Code Nr. 4512-201-0476. overall weight 4,8 kg. Materials see Recycling Passport.