





Product name	BDAS RACK
Identification code(s)	4522-130-8900*
Total weight [kg]	50 kg
Manufacturer	Philips Medical Systems Nederland B.V.
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	System_BDASRack_452213089000_579_v02

Recycling information

Special Attention	Items	Location
	Not applicable	

Batteries	Items, type	Location
	Not applicable	

Hazardous substances	Items	Location
	PWB's (up to 6 depending on model) + backpanel PWB	Figure 2,6
	Power Supply Unit front (PWB) + PSU backside (PWB) (2x)	Figure 3,6
	Frontplate different PWB's / PSU fr. (Ni layer on both sides)	Figure 1
	Frontplate BDAS (NiCr layer on both sides)	Figure 1
	Backplate BDAS (Zn Yellow passivated 5µm on both sides)	Figure 5
	Be-Cu contact springs on frontplate PWB + frontplate BDAS	Figure 1,2,4
	Remark: PWB = Printed Circuit Board.	

Fluids and gases	Items	Location
	Not applicable	

Name: Recycling passport BDAS RACK		12NC: 4522 130 89000		Version: 02 Doc.Date: 2007-05-11	
Author: G.Geboers	Page 1 of 3	Group No: 5790	Section: -	Status: Auth'd	

Locations as mentioned in the Passport (Information in pictures)



Frontplate PWB/PSU

P(rinted) W(iring) B(oard)

PSU front



Figure 1

Figure 2

Figure 3

Frontplate BDAS

BU-CU contact spring

Backplate

PSU backside





Figure 4

Figure 5

Figure 6

Backpanel

Name: Recycling passport BDAS RACK		12NC: 4522 130 89000	Version: 02 Doc.Date: 2007-05-11
Author: G.Geboers	Page 2 of 3	Group No: 5790	Section: -
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2007			Status: Auth'd

Recycling Passport Revision History

Revision	Description	Date
1	Initial version	2006-08-18

Recycling Passport Approval Signatures

Written By G. Geboers	Reviewed By (optional)
Approved By Geert-Jan Plattèl	Approved By

End of document.