





<b>Product name</b>	<b>15ch Spine coil</b>
Identification code(s)	1.5T: 4522-132-6436.0 3.0T: 4522-132-6437.0
Total weight [kg]	8.2 kg
<b>Manufacturer</b>	<b>Invivo, corporation. USA</b>
Address	3545 SW 47th Avenue
Zip code	Gainesville, FL 32608
Country	USA
Electronic info	RFCoils_Invivo15cSpineCoils_452213264360_579_v00.doc

### Recycling information

Special Attention	Items	Location
	None	Figure ..

Batteries	Items, type	Location
	None	Figure ..

Hazardous substances	Items	Location
	Several printed circuit boards inside the coil.	Figure 1
	Beryllium copper contact pins on connectors	Figure 1

Fluids and gases	Items	Location
	Phantom	Figure 2

Name: <b>Recycling Passport 15-ch Spine coil</b>		12NC: <b>4522-132-64360</b>	Version: <b>00</b> Doc.Date: <b>2007-11-05</b>
Author: <b>P. Ebben</b>	Pag.: <b>1 of 3</b>	Group No: <b>5790</b>	Section: <b>-</b> Status: <b>Authorized</b>

## Locations as mentioned in the Passport (Information in pictures)

<Place sketch(es) or photo(s) here, and indicate issues using icons >



Figure 1



Fluid inside. See label for details



Figure 2 (1.5T and 3.0T look similar)

Name: <b>Recycling Passport 15-ch Spine coil</b>		12NC: <b>4522-132-64360</b>	Version: <b>00</b> Doc.Date: <b>2007-11-05</b>
Author: <b>P. Ebben</b>	Pag.: <b>2 of 3</b>	Group No: <b>5790</b>	Section: <b>-</b> Status: <b>Authorized</b>

## Recycling Passport Revision History

Revision	Description	Date
00	Initial version	2007-11-16

## Recycling Passport Approval Signatures

<p><b>Written By</b></p> <p>Paul Ebben</p>	<p><b>Reviewed By (optional)</b></p>
<p><b>Approved By</b></p> <p>Geert-Jan Platel</p>	<p><b>Approved By</b></p>

<END OF DOCUMENT>

Name:	<b>Recycling Passport 15-ch Spine coil</b>	12NC:	<b>4522-132-64360</b>	Version: <b>00</b>	Doc.Date: <b>2007-11-05</b>
Author: <b>P. Ebben</b>		Pag.: <b>3 of 3</b>	Group No: <b>5790</b>	Section: <b>-</b>	Status: <b>Authorized</b>