




Product name	SENSE Torso / Cardiac coil
Identification code(s)	4535-300-9652x 1.5T Anterior 4535-300-9651x 1.5T Posterior 4535-301-0119x 3.0T Anterior 4535-301-0118x 3.0T Posterior 4535-302-1440x Phantom
Total weight [kg]	5.2 Kg (coils) 5 Kg (Phantom)
Manufacturer	Invivo, corporation. USA
Address	3545 SW 47th Avenue
Zip code	Gainesville, FL 32608
Country	USA
Electronic info	Recycling Passport 32-ch Cardiac.doc

Recycling information

Special Attention	Items	Location
	None	

Batteries	Items, type	Location
	None	

Hazardous substances	Items	Location
	Several printed circuit boards inside the coil	Figure 1
	Beryllium copper contact pins on connectors	

Fluids and gases	Items	Location
	Phantom fluid inside 16 small jars	Figure 2
	Fluid: 31.5mM/L of NiCL2-6H2O	

Locations as mentioned in the Passport (Information in pictures)




Figure 1



16X

Figure 2