





<b>Product name:</b>	<b>Filter Box #2: Shim RF Feedthrough</b>	
<b>Identification code</b>	4535 673 4618* (shim filter assy)	
<b>Total weight (in Kg)</b>	Approximately 40 Kg	
<b>Producer /Manufacturer</b>	<b>Name company:</b>	Philips Medical Systems (Cleveland), Inc
	<b>Address:</b>	595 Miner Rd., Cleveland, Ohio
	<b>Zip code:</b>	44143
	<b>Country:</b>	USA
	<b>Electronic info:</b>	<a href="http://www.medical.philips.com/main/company/sustainability/recycling/mri/">http://www.medical.philips.com/main/company/sustainability/recycling/mri/</a>

<b>Recycle info</b>	<b>Substances</b>	<b>Location</b>
<b>Hazardous</b>  <b>To be removed</b>	Lead(Pb) (Solder – PC boards)	Figure 2
 <b>Batteries</b> <b>To be removed</b>	<b>Type</b>	<b>Location</b>
<b>Special attention</b> 	<b>Item</b>	<b>Location</b>
 <b>Fluids / Gases</b>	<b>Item</b>	<b>Location</b>

\* Denotes revision level wildcard designation. (Numbers range from 1-9).

Locations as mentioned in the Passport (Pictures information)



Figure 2 7T Shim filter assembly (lower unit) RF feedthrough.