





Product name	OCD4 (OPERATOR CONSOLE MONITOR-FIMI-Q19)
Identification code(s)	4522 132 3316x which is part of OCD-set (incl.cables) 4522 132 34592
Total weight [kg]	8 kg
Manufacturer	FIMI
Address	Via Saul Banfi, 1 - 21047 Saronno
Zip code	
Country	Italy
Electronic info	Computer_OCD4_452213233160_579_v02.doc

Recycling information

Special Attention	Items	Location
	Printed Circuit Boards (produced by lead free soldering process)	Figure 1 (item2)

Batteries	Items, type	Location
	NO	Figure ..

Recycle info	Substances		Location
	Type	Quantity	
 To be removed	Cd	0	
	Hg	2,5mg max. per lamp(*)	Figure 1 (item 1)
	Pb	0	
	Cr ⁶⁺	0	
	PBB	0	
	PBDE	0	
	(*) Mercury is present in the 4 lamps of Backlight, according to RoHS limitation		

Fluids and gases	Items	Location
	NO	Figure ..

Locations as mentioned in the Passport (Pictures information)




FIG. 1

Monitor material		
Fe	620 g	
Al	590 g	
Cu	100 g	
Plastics	1100 g	
Boards ($S^2 > 10\text{cm}^2$)	98 cm^2 / 65 g	Scaler Board
	35 cm^2 / 15 g	Keyboard
	225 cm^2 / 270 g	Power Board
LCD Panel	2980 g	19"
Pedestal material		
Fe	1920 g	
Plastics	480 g	