



Recycling Passport

Product Name:		Philips Fluoroscopy System
Product Model		Essenta RC
Identification Code		G.ESENTA RC R2
Total weight (in Kg)		1857.302kg
Producer /Manufacturer	Name	Philips and Neusoft Medical Systems Co., Ltd
	Address:	Neusoft Park, Hun Nan Industrial Area, Shenyang 110179, China

Product Name:		Geometry Imaging Module
Identification Code		4541 103 20781
Total weight (in Kg)		5.89kg

Recycle info	Substances	Location
	Iron(Fe)	Fig 2
	Plumbum(Pb)	Fig 2

Material content (global)	Weight in Kg
Al	5.5
Cu	0.04
Fe	0.05
Printed Circuit Boards	0.3



Fig 1






Fig 2

Product Name:	Main cabinet
Identification Code	4541 102 58711
Total weight (in Kg)	100kg

Recycle info	Substances	Location
	Iron(Fe)	Fig 3
	Plumbum(Pb)	Fig 3

Material content (global)	Weight in Kg
Fe	50
Cu	0.2676
Printed Circuit Boards	3.03
All other material types	46.7






Fig 3

Product Name:	Core control cabinet
Identification Code	4541 103 71901
Total weight (in Kg)	45.062kg

Recycle info	Substances	Location
	Iron(Fe)	Fig 4 Fig 5
	Plumbum(Pb)	Fig 5

Material content (global)	Weight in Kg
Fe	45
Cu	0.002
Printed Circuit Board	0.06




Fig 4




Fig 5

Product Name:	Spot film device
Identification Code	4541 102 65561
Total weight (in Kg)	68kg

Recycle info	Substances	Location
	Iron(Fe)	Fig 6 Fig 7
	Plumbum(Pb)	Fig 6 Fig 7

Material content (global)	Weight in Kg
Fe	61.7
Al	1.3
Pb	3.6
Cu	0.0854
Printed Circuit Board	1.32




Fig 6




Fig 7

Product Name:	Stand
Identification Code	4541 103 32221
Total weight (in Kg)	1515kg

Recycle info	Substances	Location
	Iron(Fe)	Fig 8
	Plumbum(Pb)	Fig 8

Material content (global)	Weight in Kg
Fe	1400
Al	25
Rubber	10
Plastics	10
Fiber glass	50
Pb	20







Fig 8

Product Name:		HV Generator
Identification Code		454110704181 454110704191 454110704201
Total weight (in Kg)		123.35kg
Producer /Manufacturer	Name	Philips and Neusoft Medical Systems Co., Ltd
	Address:	Neusoft Park, Hun Nan Industrial Area, Shenyang 110179, China

Recycle info	Substances	Location
	Transformer oil	Fig 9
	Iron(Fe)	Fig 9
	Plumbum(Pb)	Fig 9

Material content (global)	Weight in Kg
Transformer oil	18.35
Fe	45
Al	20
Cu	10
Plastics	5
Printed Circuits Boards	8
All other material types	17




Fig 9