





Printed copies are uncontrolled

<b>Product name:</b>	<b>Gemini TF Big Bore PRS Cabinet</b>	
<b>Identification code(s)</b>	<b>4535-675-2847x, ASSY, COMPLETE PRS CABINET</b> <b>(x) Denotes revision level wildcard designation. (Numbers range from 1-9)</b>	
<b>Total weight (in Kg)</b>	<b>~36 Kg</b>	
<b>Producer/ Manufacturer</b>	Name company:	Philips Medical Systems (Cleveland), Inc.
	Address:	595 Miner Road, Cleveland, Ohio
	Zip code:	44143
	Country:	United States of America

<b>Recycle Info</b>	<b>Items:</b>	<b>Location (Qty.)</b>
<b>Special attention</b> 	None	
<b>Fluids / Gases</b> 	None	
<b>Batteries</b> 	<b>Type:</b> Lithium Ion Button Cell	Location (Qty.) See Figure 1 (1)
<b>To be Removed</b>		
<b>Hazardous</b> 	<b>Substances:</b>	<b>Location (Qty.)</b>
	Printed Circuit Board Assembly (PCBA) (*) – Computer	See Figure 1 (1)
	Printed Circuit Board Assembly (PCBA) (*) – Modules	See Figure 2 (2)
	Printed Circuit Board Assembly (PCBA) (*) – Power Supply	See Figure 2 (1)
	Fan & Motor Assembly	See Figure 2 (2)
<b>To be Removed</b>	* This cabinet contains RoHS compliant and non-RoHS compliant PCBAs and components. Refer to the markings on the component for recycling purposes.	

Title: TFBB PRS Recycling Passport  
 Recycling Passport Number: DMR104554 Rev: 01

Printed copies are uncontrolled



Figure 1, PRS Front




Figure 2, PRS Rear