


Printed copies are uncontrolled

Product name:	Gemini TF Big Bore PET/CT Host Cabinet	
Identification code(s)	4535-675-2522x, ASSY, PET-CT HOST CABINET (x) Denotes revision level wildcard designation. (Numbers range from 1-9)	
Total weight (in Kg)	~36 Kg	
Producer/ Manufacturer	Name company:	Philips Medical Systems (Cleveland), Inc.
	Address:	595 Miner Road, Cleveland, Ohio
	Zip code:	44143
	Country:	United States of America

Recycle Info	Items:	Location (Qty.)
Special attention 	None	
Fluids / Gases 	None	
Batteries 	Type: Lithium Ion Button Cell	See Figure 1 (1)
To be Removed		
Hazardous 	Substances:	Location (Qty.)
	Printed Circuit Board Assembly (PCBA) (*) – Computer	See Figure 1 (1)
	Printed Circuit Board Assembly (PCBA) (*) – Module	See Figure 2 (1)
	Printed Circuit Board Assembly (PCBA) (*) – Power Supply	See Figure 2 (2)
	Fan & Motor Assembly	See Figure 2 (2)
	Printed Circuit Board Assembly (PCBA) (*) – Ethernet Switch	See Figure 2 (1)
* This cabinet contains RoHS compliant and non-RoHS compliant PCBAs and components. Refer to the markings on the component for recycling purposes.		


Title: TFBB PET/CT Host Recycling Passport
 Recycling Passport Number: DMR104555 Rev: 01

Printed copies are uncontrolled



-  Lithium Ion Battery
-  Computer

Figure 1, PET/CT Host Front



-  PCBA - Ethernet Switch
-  PCBA - Power Supply
-  PCBA - Module
-  Fan & Motor

Figure 2, PET/CT Host Rear