





Recycling Passport

Product name:	Practix 360	
Identification code	9890-010-88401	
Order No.		
Device Type	B2B (Business to business)	
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH
	Address:	Roentgenstrasse 24
	Zip code:	22335 Hamburg
	Country:	Germany
	Web info:	See Philips Intranet
Total weight (in Kg) thereof		170 Kg

Klasse	Merkmalbezeichnung	Gesamtgewicht	Gewichtseinheit
steel, iron	cast iron, cast steel	103,445	KG
steel, iron		103,445	KG
nonferrous metals and alloys	aluminum, -alloy	8,6	KG
nonferrous metals and alloys		8,6	KG
glass / ceramics	ferrite	0,18	KG
glass / ceramics		0,18	KG
plastics / organic substances	thermoplastic	10	KG
	oil	7	KG
plastics / organic substances		17	KG
batteries	batteries, lithium	0,009	KG
batteries		0,009	KG
standard parts	other electronic powered devic	10,945	KG
	computers and accessories	9,4	KG
	printed circuit boards	6	KG
	wheels / castors	5,5	KG
	cables	5	KG
	X-ray tubes	1,2	KG
	power supplies / transformers	0,58	KG
	compensation filters	0,24	KG
standard parts		38,865	KG
relevant materials	confirmed by supplier		
relevant materials			
prohibited materials	lead and -compounds	2	KG
	confirmed by supplier		
prohibited materials		2	KG
		170,099	KG

Printed copies are uncontrolled documents

Recycling Passport


Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances  To be removed	Substances:	Location
	Lead	Fig. 1 (1)
Batteries  To be removed	Type:	Location
	Lithium Battery	Fig. 1 (2)
Special attention 	Item	Location
	Discharge Power Capacitors before handling the unit	Fig. 1 (3)
Fluids / Gases 	Item	Location
	Oil tank / Remove oil before dismantling	Fig. 1 (1)
	HV-transformer oil contains no PCBs	

Printed copies are uncontrolled documents

Recycling Passport

Locations as mentioned in the Passport (Pictures information)

Fig. 1



Printed copies are
uncontrolled documents