





Recycling Passport

Product name:	Wireless Portable Detector		
Identification code	9890-010-87751		GXR-DE-3685/09
Order No.			
Device Type	B2B (Business to business)		
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	See Philips Intranet	
Total weight (in Kg) thereof			28,6 Kg

Class	Material description	Weight	Unit
steel, iron	cast iron, cast steel	6,782	KG
	iron, low alloy (<5%)	0,4	KG
steel, iron		7,182	KG
nonferrous metals and alloys	aluminum, -alloy	2,631	KG
	magnesium, -alloy	1,659	KG
	molybdenum, -alloy	0,071	KG
	zinc, -alloy	0,048	KG
	copper, -alloy	0,018	KG
nonferrous metals and alloys		4,427	KG
glass / ceramics	glass	0,445	KG
glass / ceramics		0,445	KG
plastics / organic substances	thermoplastic	5,483	KG
	thermoset	0,365	KG
	elastomer	0,006	KG
plastics / organic substances		5,854	KG
batteries	BATTERIE_LITHIUM-ION	0,286	KG
batteries		0,286	KG
standard parts	cables	6,433	KG
	printed circuit boards	1,795	KG
	power supplies / transformers	0,8	KG
	rasters	0,399	KG
	compensation filters	0,299	KG
	motors, pumps	0,012	KG
	mounting parts, attaching part	0,001	KG
standard parts		9,739	KG
relevant materials	lead and -compounds (cables an	0,391	KG
	thallium and -compounds	0,001	KG
	confirmed by supplier		
relevant materials		0,392	KG
prohibited materials	lead and -compounds	0,25	KG
	confirmed by supplier		
prohibited materials		0,25	KG
		28,575	KG

Printed copies are uncontrolled documents

Recycling Passport

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances  To be removed	Substances:	Location
	Lead	Fig. 1
	thallium	Fig. 1
Batteries  To be removed	Type:	Location
Special attention 	Item	Location
Fluids / Gases 	Item	Location

Printed copies are uncontrolled documents

Recycling Passport

Locations as mentioned in the Passport (Pictures information)

Grid frame



1

Wireless portable detektor



Dockingstation



Printed copies are
uncontrolled documents