





Recycling Passport

Product name:	Practix 300/400		
Identification code	9890-010-82692		
Order No.			
Device Type	B2B (Business to business)		
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	http://www.medical.philips.com/main/company/sustainability/	
Total weight (in Kg)			284 Kg
thereof			

Klasse	Merkmalbezeichnung	Gesamtgewicht	Gewichtseinheit
steel, iron	iron, low alloy (<5%)	41,037	KG
	iron, high alloy (>5%)	0,013	KG
	steel, iron	41,05	KG
nonferrous metals and alloys	aluminum, -alloy	89,6	KG
	copper, -alloy	0,006	KG
nonferrous metals and alloys		89,606	KG
glass / ceramics	ferrite	0,1	KG
glass / ceramics		0,1	KG
plastics / organic substances	oil	19,09	KG
	thermoplastic	7,61	KG
	elastomer	0,761	KG
plastics / organic substances		27,461	KG
batteries	batteries, lithium	0,003	KG
batteries		0,003	KG
standard parts	other electronic powered devic	48,487	KG
	X-ray tubes	17,28	KG
	cables	11,146	KG
	mounting parts, attaching part	0,798	KG
	computers and accessories	0,218	KG
	printed circuit boards	0,124	KG
standard parts		78,053	KG
relevant materials	lead and -compounds (cables an	48	KG
relevant materials		48	KG
		284,273	KG

Printed copies are uncontrolled documents


Recycling Passport

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances  To be removed	Substances:	Location
	Lead	Fig. 1
Batteries  To be removed	Type:	Location
	Lithium Battery	Fig. 2
Special attention 	Item	Location
Fluids / Gases 	Item	Location
	Oil tank / Remove oil before dismantling	Fig. 1
	HV-transformer oil contains no PCBs	

Printed copies are uncontrolled documents

Recycling Passport

Locations as mentioned in the Passport (Pictures information)



Printed copies are
uncontrolled documents