





Product name	Small Extremity Coil
Identification code(s)	9896-032-02161 1.5T 9896-032-00501 3.0 T
Total weight [kg]	1.814 kg
Manufacturer	Invivo
Address	3545 SW 47 th Ave., Gainesville, FL
Zip code	32608
Country	USA
Electronic info	SEC passport.doc

Recycling information

Special Attention	Items	Location
	none	

Batteries	Items, type	Location
	none	

Hazardous substances	Items	Location
	Several circuit boards within the coil	Figure 1
	Beryllium copper contact pins on connectors	

Fluids and gases	Items	Location
	Plastic bottle containing 2 Liters of	Figure 2
	.0315 Molal Nickel Chloride	

Name: Recycling Passport 1.5/3.0T Small Extremity Coil		12NC: 9896-032-02161 9896-032-00501		Version: 01 Doc.Date: 2008-04-101	
Author:		Pag.: 1 of 2	Group No: 5790	Section: -	Status:

Locations as mentioned in the Passport (Information in pictures)



Figure 1



Figure 2

Name: Recycling Passport 1.5/3.0T Small Extremity Coil		12NC: 9896-032-02161 9896-032-00501		Version: 01 Doc.Date: 2008-04-101	
Author:		Pag.: 2 of 2	Group No: 5790	Section: -	Status: