



| | |
|------------------------|--|
| Product name | Philips 3.0T 32-Channel Digital Head Coil |
| Identification code(s) | 9896-032-06931 |
| Total weight [kg] | 10.43 kg |
| Manufacturer | Invivo |
| Address | 3545 SW 47 th Ave., Gainesville, FL |
| Zip code | 32608 |
| Country | USA |
| Electronic info | 208161.doc |

Recycling information

| Special Attention | Items | Location |
|---|-------|----------|
|  | None | N/A |
| | | |
| | | |

| Batteries | Items, type | Location |
|---|-------------|----------|
|  | None | N/A |
| | | |
| | | |

| Hazardous substances | Items | Location |
|---|---|----------|
|  | Circuit boards within the coil with lead solder | Figure 1 |
| | Beryllium copper on contact pins on connectors | |
| | and three springs in anterior latch assemblies | |

| Fluids and gases | Items | Location |
|---|-------|----------|
|  | None | N/A |
| | | |
| | | |

| | | | | | |
|---|--|--------------------------------|--|----------------------|--|
| Name: Recycling Passport Philips 3.0T 32-CH Digital Head Coil | | 12NC: 9896-302-06931 | | Doc.Date: 5/23/2012 | |
| Pag.: 1 of 2 | | Group No: | | Section: - | |
| | | | | Status: Rev 1 | |

Locations as mentioned in the Passport (Information in pictures)



Figure 1

| | | | | | |
|---|--|--------------------------------|--|----------------------|--|
| Name: Recycling Passport Philips 3.0T 32-CH Digital Head Coil | | 12NC: 9896-302-06931 | | Doc.Date: 5/23/2012 | |
| Pag.: 2 of 2 | | Group No: | | Section: - | |
| | | | | Status: Rev 1 | |