





<b>Product name</b>	<b>Invivo iRF Foot/Ankle coils</b>
Identification code(s)	9896-032-05861 1.5T 8Ch FAC 9896-032-05871 3.0T 8Ch FAC
Total weight [kg]	7 kg
<b>Manufacturer</b>	<b>Invivo, Corporation. USA.</b>
Address	3545 SW 47th Avenue
Zip code	Gainesville, FL 32608
Country	USA
Electronic info	<a href="#">RFCoils_Invivo_iRF_FAC_Coil.pdf</a>

### Recycling information

Special Attention	Items	Location
	Several lead-based soldered printed circuit boards inside the coil.	Figure 1 (1)

Batteries	Items, type	Location
	None.	

Hazardous substances	Items	Location
	Several lead-based soldered printed circuit boards inside the coil.	Figure 1 (2)
	Beryllium copper contact pins on connectors between cable and coil.	Figure 1 (3)
	Plastic bottle with Nickel Chloride 31.5mM/L of NiCL <sub>2</sub> -6H <sub>2</sub> O.	Figure 2 (1)

Fluids and gases	Items	Location
		

Name: Recycling Passport <Component name>	12NC: <b>1234 123 1234.0</b>	Version: <b>00</b> Doc.Date: <b>2008-04-10</b>
Author: <Your Name>	Pag.: <b>1 of 2</b>	Group No: 5790    Section: -    Status: <b>Draft</b>

## Locations as mentioned in the Passport (Information in pictures)




Figure 1.



Figure 2.

Name: Recycling Passport <Component name>		12NC: <b>1234 123 1234.0</b>		Version: <b>00</b> Doc.Date: <b>2008-04-10</b>	
Author: <Your Name>		Pag.: <b>2 of 2</b>	Group No: 5790	Section: -	Status: <b>Draft</b>
© PHILIPS HEALTHCARE NEDERLAND BV, BEST 2008			Ref. XJV032011.07 Rev.01		