





Product name	Invivo iRF Wrist coils
Identification code(s)	9896-032-05881 1.5T 8Ch Wrist 9896-032-05891 3.0T 8Ch Wrist
Total weight [kg]	7 kg
Manufacturer	Invivo, Corporation. USA.
Address	3545 SW 47th Avenue
Zip code	Gainesville, FL 32608
Country	USA
Electronic info	RFCoils_Invivo_iRF Wrist_Coil.pdf

Recycling information

Special Attention	Items	Location
	Several lead-based soldered printed circuit boards inside the coil.	Figures 1 (1)

Batteries	Items, type	Location
	None.	

Hazardous substances	Items	Location
	Several lead-based soldered printed circuit boards inside the coil.	Figures 1 (2)
	Plastic container with Nickel Chloride. 31.5mM/L of NiCL ₂ -6H ₂ O.	Figures 1 (3)
	Beryllium copper contact pins on connectors between cable and coil.	Figures 1 (4)

Fluids and gases	Items	Location
		

Name: Recycling Passport <Component name>	12NC: 1234 123 1234.0	Version: 00 Doc.Date: 2008-04-10
Author: <Your Name>	Pag.: 1 of 2	Status: Draft
	Group No: 5790	Section: -

Locations as mentioned in the Passport (Information in pictures)




Figure 1



1



2



3

Name: Recycling Passport <Component name>		12NC: 1234 123 1234.0	Version: 00 Doc.Date: 2008-04-10
Author: <Your Name>	Pag.: 2 of 2	Group No: 5790	Section: - Status: Draft