





<b>Product name</b>	<b>OCD4 (OPERATOR CONSOLE MONITOR190P7)</b>
Identification code(s)	4522 132 3316x which is part of OCD-set (incl.cables) 4522 132 34591
Total weight [kg]	5.67 kg
<b>Manufacturer</b>	<b>Philips CE ; Philips Nederland -PC Peripherals</b>
Address	Vredeoord 105; Postbus 90050
Zip code	5621 CX Eindhoven
Country	The Netherlands
Electronic info	Computer_OperatorsConsoleDisplay_452213234590_579_v03.doc

## Recycling information

<b>Special Attention</b>	<b>Items</b>	<b>Location</b>
	Printed Circuit Boards	Figure 1


<b>Batteries</b>	<b>Items, type</b>	<b>Location</b>
	NO	Figure ..

<b>Hazardous substances</b>	<b>Items</b>	<b>Location</b>
		Figure 1
	Mercury contained in backlight which are exempted from RoHS directive.	Figure 2

<b>Fluids and gases</b>	<b>Items</b>	<b>Location</b>
	NO	Figure ..


Name: <b>Recycling Passport &lt;Component name&gt;</b>		12NC: <b>4522 132 34590</b>		Version: <b>03</b> Doc.Date: <b>2007-07-31</b>	
Author: <b>Johnson Lu</b>	Pag.: <b>1 of 3</b>	Group No: <b>5790</b>	Section: <b>-</b>	Status: <b>Authorized</b>	

**Figure 1:**



Name: <b>Recycling Passport &lt;Component name&gt;</b>		12NC: <b>4522 132 34590</b>		Version: <b>03</b> Doc.Date: <b>2007-07-31</b>	
Author: <b>Johnson Lu</b>		Pag.: <b>2 of 3</b>	Group No: <b>5790</b>	Section: <b>-</b>	Status: <b>Authorized</b>

**Figure 2:**



Mercury contained in  
backlight hidden behind  
lcd screen

### Recycling Passport Revision History

Revision	Description	Date
1	Initial version	Aug-11'06

### Recycling Passport Approval Signatures

<b>Written By</b> Johnson Lu	<b>Reviewed By (optional)</b> D.Sinnema
<b>Approved By</b> Amos Chen	<b>Approved By</b> D.Sinnema

End of document.

Name: <b>Recycling Passport &lt;Component name&gt;</b>		12NC: <b>4522 132 34590</b>	Version: <b>03</b> Doc.Date: <b>2007-07-31</b>
Author: <b>Johnson Lu</b>	Pag.: <b>3 of 3</b>	Group No: <b>5790</b>	Section: <b>-</b>
Status: <b>Authorized</b>			Template version 2.