



Optimizing CPAP prescription pressure at home

Your Home Sleep Testing and long-term CPAP solution

With System One REMstar Pro and REMstar Auto you now have the ability to assess your home sleep testing patients CPAP prescription in an Auto-CPAP Trial mode that was designed to optimize pressure settings, while providing them the long-term benefits of fixed CPAP therapy. This therapy uniquely adjusts to your patient's CPAP pressure needs over time through the CPAP-Check mode. EncoreAnywhere gives you the ability to remotely monitor AHI, leak and many other advanced events recorded by your patient's REMstar Pro and Auto devices.

Enhancing patient compliance

CPAP-Check mode checks on the user every 30 hours to determine if therapy pressure is optimal. If not, it automatically adjusts the nightly fixed CPAP pressure by 1 cm H₂O (with maximum upper and lower limits of 3 cm) to obtain an ideal pressure.

CPAP-Check benefits:

- Therapeutic pressure is continuously monitored
- Long-term enhanced patient care – patients are always on optimal CPAP pressure
- Therapy adapts to patient's pressure needs over time
- Less unnecessary intervention thanks to remote management via EncoreAnywhere

Auto-Trial mode delivers breath-by-breath Auto-CPAP therapy for up to a total of 30 days. At the end of the Auto-Trial period, the device will automatically transition into the CPAP-Check mode at the pressure the patient was at or below 90% of the trial period time.

Auto-Trial benefits:

- Provides clinicians the opportunity to provide short term Auto-CPAP for patients who would benefit without having to prescribe a "full" Auto-CPAP device
- Reevaluation of therapeutic pressures for patients who:
 - received a home sleep test
 - may have had a poor in-lab titration study
 - who experience a return of symptoms

¹ Available in REMstar Pro and REMstar Auto devices

PHILIPS

RESPIRONICS

sense and simplicity

System One: a quality choice

In addition to the enhanced features and appearance of System One, we've built even more quality and durability while providing you the same exceptional product and sales support. When you add it all up, it's clear why System One remains an intelligent choice.

System One quality comes through

- 25% quieter than its closest competitor
- Durable, easy to connect, top-mount tubing connections for an exceptional customer experience
- Commitment to quality excellence has earned Philips Respironics the prestigious 2011 Manufacturing Plant of the Year by *Assembly Magazine*
- Flex Promise program
- One Touch warranty program

Try one of our newest masks!

Choose TrueBlue – the Philips gold standard for comfort, stability, and seal.



Ordering information	Product number
REMstar Pro C-Flex+	DS460S
REMstar Pro C-Flex+ with heated humidifier	DS460HS
REMstar Pro C-Flex+ with heated tube humidifier	DS460TS
REMstar Auto A-Flex	DS560S
REMstar Auto A-Flex with heated humidifier	DS560HS
REMstar Auto A-Flex with heated tube humidifier	DS560TS
Accessories	Part number
System One 60 Series heated humidifier	DS6H
60 watt power supply (for S and HS models only)	1091398
System One 60 Series heated tube humidifier	DS6T
80 watt power supply (for TS model only)	1091399
System One heated tube (for TS model only)	SYS1HT15
Universal water chamber (for 50 or 60 series devices)	1063785
Power cord (for 60 series devices only)	1005894

Respironics, System One, REMstar Pro, REMstar Auto, C-Flex+, Auto A-Flex, EncoreAnywhere, CPAP-Check, Auto-Trial and TrueBlue are trademarks of Koninklijke Philips Electronics N.V. and its affiliates. All rights reserved.

Please visit www.philips.com/systemone



© 2012 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

RW 4/26/12 MCI 4104897

Philips Healthcare is part of Royal Philips Electronics

Philips Respironics International Headquarters
+33 1 47 28 30 82
Philips Respironics Asia Pacific
+65 6882 5282
Philips Respironics Australia
+61 (2) 9666 4444
Philips Respironics United Kingdom
+44 800 1300 845

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)
www.philips.com/respironics