



System One
Sleep Therapy Systems

Category: CPAP

REMstar Plus



DS260S REMstar Plus C-Flex
DS260HS REMstar Plus C-Flex
with humidifier
DS260TS REMstar Plus C-Flex
with Heated Tube
humidifier

A proven CPAP sleep therapy device with C-Flex pressure relief

System One REMstar Plus offers the revolutionary comfort of C-Flex with patterns of use data reporting to help meet your patients' compliance needs. At the same time, REMstar Plus includes virtually all of the features and programs that make System One the most prescribed sleep therapy system by US sleep centers.

System One features

- C-Flex for comfort
- Oximetry module
- On-board memory
- Encore patient management system
- SD card and wired, wireless, and broadband modem capable

Humidification

- System One humidification control
- Pre-heat humidifier option
- Heated Tube humidifier option

Programs

- Flex Promise program
- Bi-level Rescue program
- System "One-Touch" program
- Partners in Compliance Management

REMstar Plus advantages

- CPAP mode
- C-Flex for comfort
- Patterns of use data reporting
- Modem and SD card capable
- Optional, integrated oximetry available
- Heated Tube option

PHILIPS
RESPIRONICS

REMstar Plus specifications

	Comfort features	Encore data capture and reporting	
Modes	Flex comfort	Patterns of use	Reimbursement
CPAP	C-Flex	Standard	E0601

Specifications	
Pressure range	4-20 cm H ₂ O
Flex pressure relief	0 to 3
Ramp time	0 to 45 min (5-min increments)
Starting ramp pressure	4 to CPAP; patient adjustable
Dimensions	7" L x 5.5" W x 4" H 18 cm L x 14 cm W x 10 cm H
Weight	3 lbs; 1.36 kg (without humidifier)
Filters	Foam and optional ultra-fine
Device set-up	LCD/control wheel/push button
Data storage capacity (minimum)	Display: compliance VIC, 1-, 7- and 30-day averages SD card > 1 year On-board: up to 6 months
Compliance meter	Breathing detection
Altitude compensation	Manual
Electrical requirements	100-240 VAC, 50/60 Hz
Humidification	Integrates with System One 60 Series heated humidifier and System One 60 Series Heated Tube humidifier
Warranty	Two years (US only)
Miscellaneous	Fine pressure calibration

Sleep therapy accessories	
Humidification	Part number
System One 60 Series humidifier	DS6H
System One 60 Series Heated Tube humidifier	DS6T
Water chamber (includes both 50 and 60 Series lids)	1063785
Heated Tube humidifier upgrade kit	6T-CONVKIT
Tubing	Part number
System One Heated Tube	SYS1HT15
System One Performance Tube (15mm)	1074400
Performance Tubing (22mm)	1032907
Filters	Part number
Pollen (foam) filter, reusable (1/pack)	1035443
Pollen (foam) filter, reusable (2/pack)	1029330
Ultra-fine filter, disposable (1/pack)	1063096
Ultra-fine filter, disposable (2/pack)	1035442
Ultra-fine filter, disposable (6/pack)	1029331
Accessories	Part number
Device side panel	1099592
Carrying case	1063857



Complies with the Buy American Act

Respironics, Auto-Trial, C-Flex, CPAP-Check, Encore, REMstar, and System One are trademarks of Koninklijke Philips Electronics N.V. and its affiliates. All rights reserved.

Please visit www.philips.us/remstarplus



© 2012 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

edoc JJ 10/30/12 MCI 4105270

Philips Healthcare is part of Royal Philips Electronics

Philips Respironics International Headquarters
+33 1 47 28 30 82
Philips Respironics Asia Pacific
+65 6882 5282
Philips Respironics Australia
+61 (2) 9947 0440
Philips Respironics United Kingdom
+44 800 1300 845

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)
www.philips.com/respironics