



Like a glove

ComfortGel nasal mask provides a high level of fit and comfort from the moment you slip it on.

This gel mask combines premium comfort features with a host of innovations. With the ComfortGel, it's easier to get the right fit the first time – even on hard-to-fit patients.

Establishing a good seal and comfortable fit is a critical part of successful long-term therapy. The ComfortGel provides both for men and women wearers. We've paired our proprietary blue gel with a proven dual-layer cushion design for an unparalleled seal. The adjustable StabilitySelector simplifies fitting and optimizes seal, while the premium headgear with convenient ball-and-socket headgear attachments makes it easy to remove the mask. And getting a good fit quickly is essential to your business. It all adds up to one streamlined mask that helps you fit 100% of your patients, 100% of the time.

PHILIPS

RESPIRONICS

Removable forehead support and pad adjust to match the angle of the wearer's forehead.

Adjustable StabilitySelector ensures best mask seal and fit.

360° swivel allows free movement of tubing.

"Silent" exhalation port means quiet operation.

Ball-and-socket headgear attachments enhance movement and allow for easy mask removal.



ComfortGel uses one-size-fits-most premium headgear with EZ Peel tabs, which adjust easily in four places for a secure, custom fit. The headgear evenly distributes tension, adds stability, and helps to prevent mask leaks. ComfortGel disassembles easily for cleaning and replacement. The headgear can be used with most of our masks.

Replaceable silicone Comfort Flap improves seal and simplifies cleaning.

Inner blue gel layer provides comfort and support and can be easily customized.

Item					
ComfortGel mask	Part number	HCPCS code	DuoPack (Extra gel and silicone cushions)		
			Part number	HCPCS code	
Petite mask with headgear	1009040	A7034,A7035	1035999	A7034,A7035,A7032	
Small mask with headgear	1009041	A7034,A7035	1036000	A7034,A7035,A7032	
Medium mask with headgear	1009042	A7034,A7035	1035931	A7034,A7035,A7032	
Large mask with headgear	1009043	A7034,A7035	1035932	A7034,A7035,A7032	
Petite mask only	1010518	A7034	N/A	N/A	
Small mask only	1010519	A7034	N/A	N/A	
Medium mask only	1010640	A7034	N/A	N/A	
Large mask only	1010641	A7034	N/A	N/A	
FitPack: medium mask with headgear plus small gel and silicone cushions	1031405	A7034,A7035	N/A	N/A	
FitPack: medium mask without headgear plus small gel and silicone cushions	1031406	A7034	N/A	N/A	
Replacement cushions	Petite	Small	Medium	Large	HCPCS code
Gel cushion	1009048	1009049	1009050	1009051	A7032
Silicone Comfort Flap	1009052	1009053	1009054	1009109	N/A
ComfortGel cushion and silicone Comfort Flap	1031391	1031392	1031403	1031404	A7032
Parts and accessories	Part number	HCPCS code	Parts and accessories	Part number	HCPCS code
Premium headgear with EZ Peel tabs	1033678	A7035	Exhalation elbow	1010646	A7045
Forehead support and pad	1018401	N/A	Mask swivel	7041	N/A
Silicone forehead pads (4 per pack)	1009722	N/A	Sizing gauge	1012040	N/A
Ball-and-socket Quick Clips (2 per pack)	1009723	N/A			

Respironics, Comfort Flap, ComfortGel, DuoPack, EZ Peel, FitPack, Quick Clip and StabilitySelector are trademarks of Respironics, Inc. and its affiliates.



© 2010 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Hoeh JJ 3/2/10 MCI 4103299 PN 1044392

Philips Healthcare is part of Royal Philips Electronics

Respironics Asia Pacific
+65 6882 5282
Respironics Australia
+61 (2) 9666 4444
Respironics Europe, Middle East, Africa
+33 1 47 52 30 00
Respironics United Kingdom
+44 800 1300 845

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)
www.philips.com/respironics