

Lumify Device Compatibility List - China



Supported Devices- China 兼容设备列表*

Device Name -型号 / 号	Screen Size 尺寸	Connectivity Options 无线方式	Internal Storage 内存	Approximate Continuous Scan Time - 连续扫查时间	USB Type 探头接口
Huawei Mediapad M5 (Model SHT-W09, SHT-AL09)	8.4 in	Wi-Fi or Wi-Fi/LTE	64 GB	3 h 30 min	C
Nokia 6 (Model TA-1003)	5.5 in	Wi-Fi/LTE	32 GB	2 h 30 min	B
Huawei Mediapad M3 (Model BTV-W09, BTV-DL09)	8.4 in	Wi-Fi or Wi-Fi/LTE	32 GB	3 h 30 min	B
Samsung Galaxy Tab S2 9.7" (Only Model SM-T813)	9.7 in	Wi-Fi	32 GB	3 h	B
Samsung Galaxy Tab S2 8" (Only Model SM-T713)	8 in	Wi-Fi	32 GB	2 h 30 min	B
Panasonic Toughpad FZ-A2	10.1 in	Wi-Fi/LTE	32 GB	3 h 30 min	C
Panasonic Toughpad FZ-B2 (USB C models only)	7 in	Wi-Fi/LTE	32 GB	3 h 30 min	C
Samsung Galaxy Tab S4 (Model SM-T830)	10.5 in	Wi-Fi	64 GB	5 h 30 min	C
Samsung J7 (Model SM-J710)	5.5 in	Wi-Fi/LTE	16 GB	4 h 15 min	B
Google Pixel 2 XL	6.0 in	Wi-Fi/LTE	64 GB	3 h	C
Samsung Galaxy S9 (SM-G960U1/ SM-G9600)	5.8 in	Wi-Fi / LTE	32 GB	3 h 45 min	C
Samsung Galaxy S9+ (SM-G965U1)	6.2 in	Wi-Fi / LTE	64 GB	4 h	C
Samsung Galaxy S7 Edge (SM-G935)	5.5 in	Wi-Fi/LTE	32 GB	2 h 15 min	B
Samsung Galaxy Note8 (Model SM-N950U1 / SM-N9500)	6.3 in	Wi-Fi / LTE	64 GB	3 h 30 min	C
Samsung Galaxy S8 (Model SM-G950U1 / SM-G950W)	5.8 in	Wi-Fi / LTE	32 GB	3 h 30 min	C
Samsung Galaxy S8+ (Model SM-G955U1 / SM-G955W)	6.2 in	Wi-Fi / LTE	64 GB	4 h	C



Philips Ultrasound, Inc.

22100 Bothell Everett Hwy

Bothell, WA 98021-8431, USA

www.healthcare.philips.com/ultrasound

© Koninklijke Philips N.V. 2017. All rights reserved. Printed in USA.

Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical, or otherwise, is prohibited without the prior written consent of the copyright owner.

261704 Rev A April 2019

Lumify Device Compatibility List - China



*Devices on the supported list satisfy all of the following criteria:

1. Passed software testing on the latest Lumify software,
2. Passed Philips hardware testing,
3. Is readily available from the device manufacturer. Philips cannot test every device available on the market, so this list represents a range of devices to satisfy a variety of user needs with respect to form factor, price, battery life and connectivity options



Philips Ultrasound, Inc.

22100 Bothell Everett Hwy

Bothell, WA 98021-8431, USA

www.healthcare.philips.com/ultrasound

© Koninklijke Philips N.V. 2017. All rights reserved. Printed in USA.

Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical, or otherwise, is prohibited without the prior written consent of the copyright owner.

261704 Rev A April 2019