

SureSigns VM4

Just what you need in a portable bedside monitor



SureSigns VM4.

Your days are already very busy. You are doing more with less, and that's not likely to change anytime soon. Trust Philips to give you what you need to help care for your patients. The SureSigns VM4 has the complete solution: a monitor with ECG and vital signs measurements in one easy-to-use, compact package.

Every Philips SureSigns monitor blends vital measurement and monitoring capabilities with convenience and reliability. Basic patient monitoring—and quality patient care—have never been more clear.





The Philips VSV allows you to check the status of up to 12 patients at a glance.

Easy to learn and easy to love

The VM4, like all SureSigns monitors, is easy to set up, easy to learn, easy to use, easy to support, and, above all, easy to love. Measurements include ECG, pulse oximetry, non-invasive blood pressure, and predictive temperature (optional). Intuitive navigation puts your most commonly used features at your fingertips. The crisp color display is bright and easy to read, even in low-light situations. SureSigns VM4 is lightweight, yet durable, with a comfortable handle for easy mobility, and a battery life of 4 hours.

Connect to the SureSigns VSV: a simple way to be everywhere at once

Having a clear picture of each patient's status, at a single glance, can help relieve stress on care teams. The Philips SureSigns Vital Signs Viewer is an affordable, dependable, and simple way to add the benefits of central viewing and alarming without the complexity or expense of a traditional PC-based central station. Connect up to 12 SureSigns VM monitors at a time. It's one more way Philips is making it easier to free yourself to do exactly what you need to do to care for your patients.

All the support you need

The SureSigns VM4 is backed by high-quality technology and support. Software upgrades, via the USB port, are remarkably quick and easy. The sleek modular design simplifies set up and repair. Our warranty adds to your peace of mind. And all Philips monitors come with 24-hour telephone support from our global response centers*.



Doesn't miss a beat.



and Masimo.**

^{**} Adapter cable required.

Product speci	fications			
Width, height, depth Weight	W: 26 cm (10.2"), H: 22 cm (8.6"), D: 14.5 cm (5.7") 3.1 kg (6.9 lbs), excluding optional			
Screen	recorder and predictive temperature 8.4" (20.3 cm) SVGA TFT-AM LCD			
NBP	Oscillometric using stepwise deflation pressure Adult/pediatric measurement range - Systolic: 30 - 255 mmHg - Diastolic: 15 - 220 mmHg - Mean: 20 - 235 mmHg Neonatal measurement range - Systolic: 30 - 135 mmHg - Diastolic: 15 - 110 mmHg - Mean: 20 - 125 mmHg • Blood pressure accuracy: ±5 mmHg • Pulse rate accuracy: ±2 bpm or ±2% (whichever is greater)			
SpO ₂	 Meets EN ISO 9919 standard Measurement range: 0 – 100% Accuracy depends on sensor 			
ECG	3-lead; heart rate range: 15 – 300 bpm			
Temperature	 Mode: predictive or monitored Method: oral, rectal, or axillary Range: 26.7 – 43.3°C (80 – 110°F) 			
Battery	 Lithium ion, 11.1 V Operating time: 4 hours with continuous monitoring and NBPs every 15 minutes 			
Data output	HL7 format, via Ethernet port			
Patients	Adult, pediatric, neonatal			

Supplies included*	To re-order		
Adult antimicrobial NBP cuff	M4555A		
NBP interconnect tubing	M1599B		
Adult finger SpO ₂ sensor	M1196A		
3-lead ECG lead set with grabber cable	989803143171 (IEC) or 989803143181 (AAMI)		
Oral temperature probe**	989803143381		
Oral temperature probe covers (one box)**	M4823A		
Recorder paper (5 rolls)**	989803136891		

contact your Philips representative
for more information.
** Depending on model ordered.

* Additional supplies available. Please

863063	SureSigns VM4
T01	Predictive Temperature (optional)
R01	Recorder (optional)

SureSigns VSV

	For more	information,	visit	www.phili	ps.con
--	----------	--------------	-------	-----------	--------

m/SureSigns

989803144001

989803144011

989803147821

A simple way to be everywhere at once

Connect with the SureSigns Vital Signs Viewer to show status from up to 12 patient monitors at a time.

Getting started

Each monitor comes with set-up guide, instructions for use, quick reference cards, user training DVD, and service manual on CD.



Supplies for the SureSigns VM4 include a wide selection of disposable and reusable SpO₂ sensors, all designed for accuracy, performance, and durability.



Sleek roll stand with integrated basket and locking casters makes it easier to take SureSigns VM4 patient monitor where it's needed. Wall mount also available. Each is designed for easy mounting and quick release.

> buy online http://shop.medical.philips.com

> > + Available only in U.S.



863067

© 2008 Koninklijke Philips Electronics N.V.

Also available Roll stand with basket

Bar code scanner kit

Wall mount

All rights are reserved. Philips Medical Systems Nederland B.V. reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Nellcor is a registered trademark of Nellcor Puritan Bennett, Inc. Masimo is a registered trademark of Masimo Corporation. Sure Temp is a registered trademark of Welch Allyn.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare healthcare@philips.com fax: +31 40 27 64 887

Printed in The Netherlands 4522 962 31541/862 * APR 2008 Philips Healthcare Global Information Center P.O. Box 1286 5602 BG Eindhoven The Netherlands