



**PHILIPS**

dStream T/R Head

MR coils

## Designed for **special head applications**

dS T/R Head coil for 1.5T and 3.0T systems with dStream

### **Key benefits**

---

- High  $B_1$  field supports excellent spectroscopy
  - Spacious design to accommodate stereotactic frames
- 

The dStream (dS) T/R (transmit/receive) coil consists of a base, sliding coil, and head support. It is designed to provide excellent spectroscopy results due to its higher  $B_1$  field. In addition, it allows imaging of patients with stereotactic frames. Open design providing good imaging homogeneity.

# dStream T/R Head coil

Coverage	30 cm
Maximum number of channels	Single transmit Single receive
Main applications	Brain Spectroscopy Patients with stereotactic frame
Coil solution type	Transmit/Receive
Data sampling	dStream digital broadband
Coil connection	T/R interface
Other features	- Mirror providing open visual field - Allows imaging of patients with stereotactic frames



© 2015 Koninklijke Philips N.V. All rights reserved.  
Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
(Royal Philips) or their respective owners.

4522 991 13101 \* SEP 2015



**How to reach us**  
Please visit [www.philips.com/healthcare](http://www.philips.com/healthcare)  
Please visit [www.philips.com/Ingenia](http://www.philips.com/Ingenia)  
[healthcare@philips.com](mailto:healthcare@philips.com)