





Recycling Passport

Product name:	DigitalDiagnost Userstation Eleava Workspot	
Identification code	9890-010-87551	
Order No.		
Device Type	B2B (Business to business)	
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH
	Address:	Roentgenstrasse 24
	Zip code:	22335 Hamburg
	Country:	Germany
	Web info:	See Philips Intranet
Total weight (in Kg) thereof		34 Kg

Version: DigitalDiagnost Userstation			
Klasse	Merkmalbezeichnung	Gesamtgewicht	Gewichtseinheit
nonferrous metals and alloys	aluminum, -alloy	1,542	KG
nonferrous metals and alloys		1,542	KG
plastics / organic substances	thermoplastic	0,018	KG
plastics / organic substances		0,018	KG
batteries	batteries, lithium	0,003	KG
batteries		0,003	KG
standard parts	computers and accessories	23,748	KG
	monitors	7	KG
	cables	0,729	KG
	printed circuit boards	0,335	KG
	power supplies / transformers	0,3	KG
	other electronic powered devic	0,2	KG
standard parts		32,312	KG
Not relevant for substances list	No relevant for subst. list		
Not relevant for substances list			
		33,875	KG

Printed copies are uncontrolled documents

Recycling Passport

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances  To be removed	Substances:	Location
Batteries  To be removed	Type:	Location
	batteries, lithium	Fig. 1
Special attention 	Item	Location
Fluids / Gases 	Item	Location
	HV-transformer oil contains no PCBs	

Printed copies are uncontrolled documents

Recycling Passport

Locations as mentioned in the Passport (Pictures information)

Userstation for DigitalDiagnost Sub-System VM ; VS; TH2



Userstation tower is not on the picture

1



Printed copies are
uncontrolled documents