

# CoughAssist T70

## Prescription

Patient name \_\_\_\_\_ Date of birth \_\_\_\_\_ Date \_\_\_\_\_

### Settings

Mode:  manual  auto

Cough-Trak:  on  off

Inspiratory pressure (range) \_\_\_\_\_ cm H<sub>2</sub>O

Expiratory pressure (range) \_\_\_\_\_ cm H<sub>2</sub>O

Frequency of treatment \_\_\_\_\_

Interface:  mask

mouthpiece

trach adaptor

Inspiratory time \_\_\_\_\_ secs.

Expiratory time \_\_\_\_\_ secs.

Pulse oximetry SpO<sub>2</sub> ≤ \_\_\_\_\_ %  NA

Titrate to achieve inspiratory pressures of 35-40 cm H<sub>2</sub>O and expiratory pressures of -35-45 cm H<sub>2</sub>O or to achieve an effective cough.

### Additional orders

### Duration of use

lifetime  other \_\_\_\_\_

### Physician information

Name (please print) \_\_\_\_\_

Signature \_\_\_\_\_

Telephone \_\_\_\_\_