

CoughAssist T70

Patient/caregiver competency checklist

Patient name: _____ Caregiver name: _____

Goals	Date goals met	Trainer
Patient/caregiver will be able to:		
1. Understand basic indications for CoughAssist T70 usage		
2. Explain the buttons, indicators, and/or connectors on the front, rear, and side panels of the device Front panel <ul style="list-style-type: none"><input type="checkbox"/> Start/stop button<input type="checkbox"/> Left button<input type="checkbox"/> Right button<input type="checkbox"/> Center button<input type="checkbox"/> Manual switch Rear and side panels <ul style="list-style-type: none"><input type="checkbox"/> AC power inlet<input type="checkbox"/> Fan exhaust<input type="checkbox"/> Remote control connector (if applicable)<input type="checkbox"/> SpO₂ connector (if applicable)<input type="checkbox"/> DC power inlet (if applicable)<input type="checkbox"/> Filter<input type="checkbox"/> Detachable battery slot (if applicable)<input type="checkbox"/> Tubing retainer<input type="checkbox"/> SD card slot<input type="checkbox"/> Breathing circuit connection		
3. Connect to AC power		
4. Assemble and connect the <ul style="list-style-type: none"><input type="checkbox"/> Patient interface<input type="checkbox"/> Tubing<input type="checkbox"/> Bacteria filter		

PHILIPS

RESPIRONICS

sense and simplicity

Goals	Date goals met	Trainer		
Patient/caregiver will be able to:				
5. Change therapy settings in the Setting access screen (if applicable) <ul style="list-style-type: none"> <input type="checkbox"/> Preset <input type="checkbox"/> Mode <input type="checkbox"/> Cough-Trak <input type="checkbox"/> Inspiratory pressure <input type="checkbox"/> Expiratory pressure <input type="checkbox"/> Inspiratory time <input type="checkbox"/> Expiratory time <input type="checkbox"/> Pause time <input type="checkbox"/> Flow 				
6. Connect to DC power and explain power source indicators <ul style="list-style-type: none"> <input type="checkbox"/> External battery <input type="checkbox"/> Detachable battery pack 				
7. Understand the information on the monitoring screen <table border="0" style="width: 100%; margin-top: 10px;"> <tr> <td style="width: 50%; vertical-align: top;"> Detailed view off <ul style="list-style-type: none"> <input type="checkbox"/> Status panel <input type="checkbox"/> Manometer panel <input type="checkbox"/> Parameters panel <input type="checkbox"/> Soft keys </td> <td style="width: 50%; vertical-align: top;"> Detailed view on <ul style="list-style-type: none"> <input type="checkbox"/> Status panel <input type="checkbox"/> Manometer panel <input type="checkbox"/> Parameters panel <input type="checkbox"/> Soft keys </td> </tr> </table>	Detailed view off <ul style="list-style-type: none"> <input type="checkbox"/> Status panel <input type="checkbox"/> Manometer panel <input type="checkbox"/> Parameters panel <input type="checkbox"/> Soft keys 	Detailed view on <ul style="list-style-type: none"> <input type="checkbox"/> Status panel <input type="checkbox"/> Manometer panel <input type="checkbox"/> Parameters panel <input type="checkbox"/> Soft keys 		
Detailed view off <ul style="list-style-type: none"> <input type="checkbox"/> Status panel <input type="checkbox"/> Manometer panel <input type="checkbox"/> Parameters panel <input type="checkbox"/> Soft keys 	Detailed view on <ul style="list-style-type: none"> <input type="checkbox"/> Status panel <input type="checkbox"/> Manometer panel <input type="checkbox"/> Parameters panel <input type="checkbox"/> Soft keys 			
8. Describe a typical treatment				
9. Demonstrate how to turn therapy on				
10. Review care, cleaning, and maintenance				

Caregiver/patient signature

Trainer signature

Respironics, CoughAssist, and Cough-Trak are trademarks of Koninklijke Philips Electronics N.V. and its affiliates. All rights reserved.



© 2012 Koninklijke Philips Electronics N.V. All rights are reserved. Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication. CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

JH 02/21/12

Philips Healthcare is part of Royal Philips Electronics

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Philips Respironics International
Headquarters
+33 1 47 28 30 82

Customer Service
+1 724 387 4000
800 345 6443 (toll free)
www.philips.com/respironics