





PHILIPS HEALTHCARE	Product Recycling Passport	Page 1 of 5
-------------------------------	-----------------------------------	-------------

Product name:	Brilliance iCT Patient Support	
Identification code(s)	4535 674 80181	
Total weight (in Kg)	395	
Producer/ Manufacturer	Name company:	Philips Medical Systems
	Address:	595 Miner Road Highland Heights, Ohio
	Zip code:	44143
	Country:	United States of America

Printed copies are uncontrolled

Recycle Info	Items:	Location
Special attention 	Removal of Motor: Couch will drop when motor is removed.	Figure 1 (1)
	Green Safety Block: When removed, couch may collapse.	Figure 1 (3)
	Giant Spring: Potential energy hazard.	Figure 1 (4)
Fluids / Gases 		
Batteries 		
To be Removed		
Hazardous 	Substances:	Location
	Printed Circuit Boards	Figure 1 & 2 (2)
To be Removed		

Show locations of materials mentioned on the previous sheet. (Picture information)

Title: Recycle Passport Brilliance iCT Patient Support
 Recycling Passport Number: 4535 675 34041 Rev: B

Printed copies are uncontrolled



Figure 1: Right Side Couch



Printed copies are uncontrolled



Figure 2: Left Side of Couch



Printed copies are uncontrolled




Figure 3: Full view of Patient Support

Title: Recycle Passport Brilliance iCT Patient Support
Recycling Passport Number: 4535 675 34041 Rev: B

PHILIPS HEALTHCARE	Product Recycling Passport	Page 5 of 5
-------------------------------	-----------------------------------	-------------

Recycling Passport Revision History

Rev	Description	Date
A	CO 2019848-Release Status A	24 February 2009
B	CO 2020759 - Removed superfluous header, added whole patient support picture and corrected Rev A description to be "CO 2019848"	09 June 2009

Printed copies are uncontrolled