






Philips Medical Systems / Product Recycling Passport

Product name:	15in TFT Touch XGA Display
Identification code	M8031A
Total weight (in Kg)	4.8 kg / 3.8 kg (pedestal: optional)
Producer /Manufacturer	Name company: Philips Medizin Systeme Boeblingen
	Address: Hewlett-Packard Strasse 2
	Zip code: Boeblingen 71034
	Country: Germany

Recycle info	Substances		Location
Hazardous	Type	Quantity	
 <b style="color: red;">To be removed	Cd	0	
	Hg	Mercury is present in the 4 lamps of backlight	Figure 1 (1.1)
	Pb	Lead is present in the soldering process of PCBs	Figure 1 (2)
	Cr ⁶⁺	0	
	PBB	0	
	PBDE	0	
Batteries	Type	Quantity	Location
 <b style="color: red;">To be removed	-	0	
Special attention	Item	Location	
	-		
Fluids / Gases	Item	Location	
	-		

Philips Medical Systems / Product Recycling Passport

Locations as mentioned in the Passport (Pictures information)



Material		
Fe	1.5 kg / 3.4 kg (*)	(*) Pedestal: optional
Al	0.1 kg	-
Cu	0.15 kg	Cables
Plastics	0.8 kg / 0.4 kg (*)	(*) Pedestal: optional
Boards (S² > 10cm²)	100 cm ² 26 cm ² 44 cm ² 35 cm ² 132 cm ² 97 cm ²	Control board (item 2.1 in the picture) Keyboard (item 2.2 in the picture) Touch-screen board (item 2.3 in the picture) Audio board (item 2.4 in the picture) Inverter board (item 2.5 in the picture) LCD driver board (item 2.6 in the picture)
LCD panel	15" / 1.4kg	Colour