





Product Recycling Passport

Product name:	4500MRI Pulse Oximeter	
Identification code(s):	3109-1 & 3109-3	
Total weight (in kg):	9 kg	
Producer /Manufacturer	Name company:	Invivo Corporation
	Address:	12601 Research Parkway Orlando, FL
	Zip code:	32826
	Country:	USA
	Electronic info:	Notice, Recycling Passport, 3109-1 & 3109-3.doc

Recycle info	Item:	Location
Special Attention	None	
		

Fluids / Gases	Type:	Location
	Not applicable	

Batteries	Substances:	Location
	Battery, Alkaline-Manganese Dioxide, 1.5V, Qty 3 (Model 3109-1 only)	Figure 2 – Item 1
	Battery, Lead-Acid, 12 V, Qty 1 (Model 3109-3 only)	Figure 3 – Item 1
To be removed		

Hazardous	Item:	Location
	Lead (Pb) – Printed circuit board solder	Figure 1 – Item 1 Figure 4 – Item 1
	Lead (Pb) – Patient Cables	Figure 5 – Item 1 Figure 6 – Item 1
	Hexavalent Chromium (Cr ⁶⁺) – Metal enclosures	Figure 1 – Item 2 Figure 3 – Item 2
	Cadmium (Cd) – Patient cables	Figure 5 – Item 1 Figure 6 – Item 1

Product Recycling Passport

Locations as mentioned in the Passport (Pictures information)



Figure 1: 4500MRI Pulse Oximeter, Models 3109-1 & 3109-3



Figure 2: 4500MRI Pulse Oximeter, Model 3109-1, Batteries

Product Recycling Passport




Figure 3: 4500MRI Pulse Oximeter, Model 3109-3, Battery



Figure 4: Power Supply (Invivo Part No. HE153A & HE153B)



Figure 5: Accessory, SpO2 Grip Sensor (Invivo Part No. 9399A)

Product Recycling Passport



Figure 6: Accessory, Oxclip MRI Finger Sensor (Invivo Part No. 9367)