





Philips Medical Systems/Product Recycling Passport

Printed copies are uncontrolled

Product name:	M3171 Access Point Controller	
Identification code(s)	453564009481	
Total weight (in Kg)	1.9	
Producer/ Manufacturer	Name company:	Philips Medical Systems
	Address:	3000 Minuteman Road, Andover, MA
	Zip code:	01810
	Country:	USA

Recycle Info	Items:	Location
 Special attention	Not Applicable	N/A
 Fluids / Gases	Items:	Location
	Not Applicable	N/A
 Batteries	Type:	Location
	Not Applicable	N/A
To be Removed		
 Hazardous	Substances:	Location
	Printed circuit board	See Figure #1
	Power Supply	See Figure #2
To be Removed		

Title: Access Point Controller Recycling Passport
 Recycling Passport Number: M3171-90333 Rev: A

Philips Medical Systems/Product Recycling Passport

Printed copies are uncontrolled




Figure #1




Figure #2