





Philip Medical systems/ Recycling Passport

Printed copies are uncontrolled

Product name:	Access Point	
Identification code(s)	453563495091	
Total weight (in Kg)	0.8 Kg	
Producer/ Manufacturer	Name company:	Philips Medical Systems
	Address:	3000 Minuteman Road, Andover.MA
	Zip code:	01810
	Country:	USA

Recycle Info	Items:	Location
Special attention 	Not applicable	
Fluids / Gases 	Not Applicable	
Batteries 	Not Applicable	
To be Removed		
Hazardous 	Substances:	Location
	Access Point Mother PC Board	See Figure 1
	Access Point Radio PC Board	See Figure 2
To be Removed		

Title: Access Point Recycling Passport
 Recycling Passport Number: _453564044701_____ Rev: _A_____

Philip Medical systems/ Recycling Passport

Printed copies are uncontrolled

Figure 1 Motherboard PC Board





Figure 2. Radio PC Board (attached to motherboard)



Title: Access Point Recycling Passport
Recycling Passport Number: _453564044701_____ Rev: _A_____