





## Philips Medical Systems/Product Recycling Passport

Printed copies are uncontrolled

<b>Product name:</b>	<b>Battery Reconditioner</b>	
<b>Identification code(s)</b>	<b>AGI-3002, M4791A</b>	
<b>Total weight (in Kg)</b>	<b>.47</b>	
<b>Producer/ Manufacturer</b>	Name company:	Philips Medical Systems
	Address:	3000 Minuteman Road, Andover, MA
	Zip code:	01810
	Country:	USA

Recycle Info	Items:	Location
<b>Special attention</b> 	Ensure that the power is removed by disconnecting the power converter from the mains supply prior to disassembly.	Figure 1
<b>Fluids / Gases</b> 	Not Applicable	
<b>Batteries</b> 	Not Applicable	
<b>Hazardous</b> 	Printed circuit board, lead solder (Pb)	Figure 2
		Figure 3

# Philips Medical Systems/Product Recycling Passport

Printed copies are uncontrolled




Figure 1




Figure 2




Figure 3

Printed circuit boards