





Philips Medical Systems/Product Recycling Passport

Product name:	CO2 Transducer Sensor	
Identification code(s)	M1460A	
Total weight (in Kg)	.17	
Producer/ Manufacturer	Name company:	Philips Medical Systems
	Address:	3000 Minuteman Road, Andover, MA
	Zip code:	01810
	Country:	USA

Recycle Info	Items:	Location
Special attention 	Internal components of the M1460A CO2 Transducer contain Beryllium Copper. Return M1460A Transducers to your local Philips Medical Systems service representative for disposal.	Figure 1
Fluids / Gases 	Items: Not Applicable	Location
Batteries 	Type: Not Applicable	Location
Hazardous  To be Removed	Substances: Not Applicable	Location

Printed copies are uncontrolled

Philips Medical Systems/Product Recycling Passport



Figure 1

Printed copies are uncontrolled