





Product (Model IP5) Recycling Passport

Product name:	Expression IP5	
Reference Identifier:	Model IP5 / REF 865471	
Total weight (in Kg):	8.43 Kg	
Producer /Manufacturer	Name company:	Invivo, a division of Philips Medical Systems
	Address:	12151 Research Parkway
	Zip code:	32826
	Country:	United States

Recycle info	Substances:	Location
Special attention 	Touch Screen Display	Figure 3 (item 4)
Fluids / Gases 	Items:	Location
Batteries 	Type:	Location
	Newsun 3v Lithium Ion battery (1x)	Figure 1 (item 1) Figure 2 (item 1)
To be removed		
Hazardous 	Item	Location
	Lead Solder	Figure 3 (item 2)
To be removed		

Additional Information as needed

Material content (global)	Weight in Kg	Materials content (continue)	Weight in Kg
Magnesium	0	Battery	0.003
Aluminum	0.179	Steel	5.89
Plastics etc	0.879	All other material types	0.677
Printed Circuits Boards	0.802		

Locations as mentioned in the Passport (Pictures information)




Figure 1




Figure 2




Figure 3



Figure 4

