





Product name:	Expression MRI Patient Monitoring System	
Identification code(s):	865214	
Total Weight (in kg):	46	
Producer / Manufacturer:	Name company:	Invivo Corporation
	Address:	12501 Research Parkway, Orlando, FL
	Zip code:	32826
	Country:	USA
	Electronic info:	A-Q2975-00153

Recycle info	Item	Location
Special Attention 	None	Not applicable
Fluids / Gases 	None	Not applicable
Batteries  To be removed	Substances Batteries, Lithium-polymer, Quantity 6 (1 battery per wireless module, 2 for the DCU, 2 for the Cart/PMC) (Invivo part numbers 9065 and 989803169491)	Location Figure 1 – Item 1 Figure 4 – Item 1 Figure 5 Figure 6 – Item 1
Hazardous  To be removed	Item Lead (Pb) – Printed circuit board solder	Location Figure 2 – Item 1 Figure 3 – Item 1 Figure 6 – Item 2

Title: Expression MRI Patient Monitoring System, Recycling Passport
 Recycling Passport Number: A-Q2975-00153 Rev: A




Figure 1: Battery Location, Wireless Modules




Figure 2: Wireless ECG Module




Figure 3: Wireless Processing Unit




Figure 4: Patient Management Configuration, Battery Removal



Figure 5: Cart Configuration, Battery Removal




Figure 6: Display Control Unit

Recycling Passport Revision History

Rev	Description	Date
A	Initial Release	4-Feb-2010

Recycling Passport Approval Signatures

Written By Krystal Mitchell, Quality & Regulatory Engineer

Reviewed By (optional) N/A

Approved By <i>Krystal Mitchell</i>

Approved By N/A
