





Philips Medical Systems/Product Recycling Passport

Printed copies are uncontrolled

Product name:	ICG Monitor	
Identification code(s)	862146, 453564012601	
Total weight (in Kg)	2.2	
Producer/ Manufacturer	Name company:	Philips Medical Systems
	Address:	3000 Minuteman Road, Andover, MA
	Zip code:	01810
	Country:	USA

Recycle Info	Items:	Location
Special attention 	Ensure that the power is removed by disconnecting the battery eliminator power cord from the mains supply prior to disassembly.	Figure 1
Fluids / Gases 	Items:	Location
	Not Applicable	
Batteries 	Type:	Location
	Sealed lead-acid (Pb)	Figure 2
Hazardous 	Substances:	Location
	Printed circuit board, lead solder (Pb)	Figure 3
	LCD backlight display (Hg)	Figure 4
To be Removed		

Philips Medical Systems/Product Recycling Passport

Printed copies are uncontrolled




Figure 1



Battery

Figure 2





Figure 3



LCD Display

Figure 4