




Philips Medical Systems / Product Recycling Passport

Product name:	MP20/MP30 Patient Monitor	
Identification code	M8001A/M8002A	
Total weight (in Kg)	5.8 Kg	
Producer /Manufacturer	Name company:	Philips Medizin Systeme Boeblingen
	Address:	Hewlett-Packard Strasse 2
	Zip code:	Boeblingen 71034
	Country:	Germany

Recycle info	Substances:	Location
Hazardous	N/A	
 To be removed		

Batteries	Type:	Location
 To be removed	Lithium Ion battery (optional: 2x, 1x, 0x)	Figure 3 (item 1)

Special attention	Item	Location
	Magnesium casting housing.	Figure 4 (item 2)
	TFT-Display	Figure 1 (item 3)

	Item	Location

Material content (global)	Weight in Kg	Materials content (continue)	Weight in Kg
Magnesium	2.5	Battery	~1
Aluminium	0.3	Steel	0.2
Plastics etc	0.9	All other material types	~0.3
Printed Circuits Boards	~0.6		

Locations as mentioned in the Passport (Pictures information)



Figure 1



3



Figure 2



Figure 3



1



Figure 4



2




Figure 5