




Philips Medical Systems / Product Recycling Passport

Product name:	MP60/MP70 Patient Monitor	
Identification code	M8005A/M8007A	
Total weight (in Kg)	10.8 Kg	
Producer /Manufacturer	Name company:	Philips Medizin Systeme Boeblingen
	Address:	Hewlett-Packard Strasse 2
	Zip code:	Boeblingen 71034
	Country:	Germany

Recycle info	Substances:	Location
Hazardous	N/A	
 To be removed		

Batteries	Type:	Location
 To be removed	N/A	

Special attention	Item	Location
	TFT-Display	Figure 1 (1)

	Item	Location

Material content (global)	Weight in Kg	Materials content (continue)	Weight in Kg
Steel	3.8		
Aluminium	1.4		
Plastics etc	2.9	All other material types	~1.1
Printed Circuits Boards	~1.5		

Locations as mentioned in the Passport (Pictures information)



Figure 1



Figure 2