

Product Recycling Passport

Product name:	Millennia Patient Monitor	
Identification code(s):	3155MVS	
Total weight (in kg):	14 kg	
Producer /Manufacturer	Name company:	Invivo Corporation
	Address:	12601 Research Parkway Orlando, FL
	Zip code:	32826
	Country:	USA
	Electronic info:	Notice, Recycling Passport, 3155MVS.doc

Recycle info	Item:	Location
Special Attention 	None	
Fluids / Gases 	Type:	Location
	Not applicable	
Batteries 	Substances:	Location
	Battery, Lithium, 3V, Qty 1	Figure 2 – Item 1
	Battery, Lead-acid, 12V, Qty 1	Figure 1 – Item 4
To be removed		
Hazardous 	Item:	Location
	Lead (Pb) – Printed circuit board solder	Figure 1 – Item 1 Figure 3 – Item 1
	Lead (Pb) – LCD Display	Figure 1 – Item 3
	Hexavalent Chromium (Cr ⁶⁺) – Metal enclosures	Figure 1 – Item 2 Figure 3 – Item 2
	Mercury (Hg) – Display Backlight	Figure 1 – Item 5
To be removed		

Product Recycling Passport

Locations as mentioned in the Passport (Pictures information)



Figure 1: 3155MVS Millennia Remote Patient Monitor



**Figure 2: 3155MVS, Visual of Power Supply Board
(located at the base of the device)**

Product Recycling Passport



Figure 3: Power Supply (Invivo Part no. AS201)