





## Philips Medical Systems/Product Recycling Passport

Printed copies are uncontrolled

<b>Product name:</b>	M3166A Power Distribution Module	
<b>Identification code(s)</b>	<b>M3166-60000, M3166A</b>	
<b>Total weight (in Kg)</b>	<b>22.23</b>	
<b>Producer/ Manufacturer</b>	Name company:	Philips Medical Systems
	Address:	3000 MinuteMan Road, Andover, MA
	Zip code:	01810
	Country:	USA

Recycle Info	Items:	Location
<b>Special attention</b> 	Remove top via 6 Phillips Head Screws	Figure 1
<b>Fluids / Gases</b> 	Items:	Location
	Not Applicable	
<b>Batteries</b> 	Type:	Location
	Not Applicable	
<b>To be Removed</b>		
<b>Hazardous</b> 	Substances:	Location
	Printed circuit board	Figure 2
<b>To be Removed</b>		

## Philips Medical Systems/Product Recycling Passport

Printed copies are uncontrolled



Figure 1 (picture of unit with cover off – removed with Phillips Screwdriver)




Figure 2 (Printed Circuit Card on one end of assembly)

Title: M3166A Power Distribution Module Recycling Passport  
Recycling Passport Number: A-M3166-90036 Rev: A