





Philips Medical Systems/ Product Recycling Passport

Printed copies are uncontrolled


Product name:	Sync Unit	
Identification code(s)	453564025901	
Total weight (in Kg)	3.4 Kg	
Producer/ Manufacturer	Name company:	Philips Medical Systems
	Address:	3000 Minuteman Road, Andover, MA
	Zip code:	01810
	Country:	USA

Recycle Info	Items:	Location
Special attention 	Not applicable	
Fluids / Gases 	Not Applicable	
Batteries 	Not Applicable	
To be Removed		
Hazardous 	Substances:	Location
	Sync PC Board	See Figure 1
	Power Supply PC Board	See Figure 1
To be Removed		

Title: Sync Unit Recycling Passport
 Recycling Passport Number: 453564044711 Rev: A

Philips Medical Systems/ Product Recycling Passport

Figure 1 Sync PC Boards



Printed copies are uncontrolled