


Philips Medical Systems/Product Recycling Passport

Printed copies are uncontrolled

Product name:	Philips Telemon B/C	
Identification code(s)	M2636B, M2636C	
Total weight (in Kg)	2.2	
Producer/ Manufacturer	Name company:	Philips Medical Systems
	Address:	3000 Minuteman Road, Andover, MA
	Zip code:	01810
	Country:	USA

Recycle Info	Items:	Location
Special attention 	Ensure that the power is removed by disconnecting the power converter from the mains supply prior to disassembly.	Figure 1
Fluids / Gases 	Not Applicable	
Batteries 	Type: Nickel Metal Hydride (NiMH)	Figure 3
Hazardous 	Substances: Printed circuit board, lead solder (Pb)	Figure 2
	LCD Display (Hg)	Figure 3
		Figure 4
To be Removed		

Philips Medical Systems/Product Recycling Passport

Printed copies are uncontrolled




Figure 1




Figure 2




Figure 3




Figure 4