





Philips Medical Systems/Product Recycling Passport

Printed copies are uncontrolled

Product name:	Universal Power Converter	
Identification code(s)	453563464761, M4792A	
Total weight (in Kg)	.44	
Producer/ Manufacturer	Name company:	Philips Medical Systems
	Address:	3000 Minuteman Road, Andover, MA
	Zip code:	01810
	Country:	USA

Recycle Info	Items:	Location
Special attention 	Ensure that the power is removed by disconnecting the power converter from the mains supply prior to disassembly.	Figure 1
Fluids / Gases 	Items:	Location
	Not Applicable	
Batteries  To be Removed	Type:	Location
	Not Applicable	
Hazardous  To be Removed	Substances:	Location
	Printed circuit board, lead solder (Pb)	Figure 2

Title: Universal Power Converter Recycling Passport
 Recycling Passport Number: A-453563464761-RP Rev: A

Philips Medical Systems/Product Recycling Passport




Figure 1




Figure 2