





**PHILIPS HEALTHCARE
PRODUCT RECYCLING PASSPORT**

Printed copies are uncontrolled

Product name:	HD15 Medical Diagnostic Ultrasound Equipment	
Identification code(s)	989605367641	
Total weight (in Kg)	81 kg	
Producer/ Manufacturer	Name company:	Philips Ultrasound
	Address:	22100 Bothell Everett Highway
	Zip code:	98012
	Country:	USA

Recycle Info	Items:	Location
Special attention 	Gas strut under pressure (2 locations)	Figure 1 [1]
Fluids / Gases 	Oil and Nitrogen gas: in gas strut in lift mechanism. Do not drain oil with device under pressure	Figure 1 [1]
Batteries 	Li Battery on CPU PCB	Figure 1 [2]
Hazardous  To be Removed	Lead (Pb) -- Contains lead in PCB's and solder.	Figure 1 [3]
	Mercury (Hg) – Monitor contains mercury.	Figure 1 [4]

Title: HD15 Recycling Passport
 Recycling Passport Number: 989605367640057 Rev: A

PHILIPS HEALTHCARE PRODUCT RECYCLING PASSPORT

Show locations of materials mentioned on the previous sheet. (Picture information)

Printed copies are uncontrolled

