



Speed, coverage and high resolution

Philips dStream Shoulder 8ch coil solution for Ingenuia 3.0T

Coil designed for high uniformity throughout the shoulder joint, with excellent penetration into the labrum. The coil consists of a base plate and an adjustable shoulder cup which can be raised and pivoted for comfortable positioning. dS Shoulder 8ch features DirectDigital sampling in the coil, capturing data where the MR signal is at its purest. Digital transmission without losses in the RF chain provides enhanced SNR and dS-SENSE enhanced parallel imaging performance.

- Digitization in the coil for increased SNR
- Adjustable design for a comfortable fit for either left or right shoulder
- High signal uniformity with excellent penetration into labrum

PHILIPS

dS Shoulder 8ch coil solution

Coverage	12 cm LR
Max. number of channels	8
Main applications	Shoulder
Coil solution type	Dedicated
Data sampling	DirectDigital
Coil connection	Single FlexConnect
SENSE parallel imaging	dS-SENSE enhanced parallel imaging performance



Please visit www.philips.com/MR



© 2014 Koninklijke Philips N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips

www.philips.com/healthcare
healthcare@philips.com

Printed in The Netherlands
4522 962 65861 * MAR 2014