



Powerful patient management on your desktop

EncorePro 2 patient management system

The new EncorePro 2 is a simple, powerful, efficient and secure desktop patient management system that helps you manage your patients and business more efficiently than before. Incorporating the best features of our EncorePro data management system—such as data downloading, compliance reporting and prescription management—EncorePro 2 adds even more features and efficiencies.

- User friendly and intuitive interface is simple to learn and use
- Notifications aid identification of patients requiring attention
- Reminders, notes, contacts and business reports help manage clinical and business workflow
- Highly secure system helps you comply with your requirements under HIPAA
- Supports all Respironics therapy devices



Recommended system requirements

Windows® XP Professional SP2, XP Home SP2, Vista® Business SP1, Vista Home SP1, Vista Ultimate SP1, or Windows Server® 2003 SP2

Microsoft .NET framework 3.5 SP1 or later

Windows Installer 4.5 or later

Microsoft® SQL Server® 2008

1 GHZ or faster processor speed

1 GB of RAM

40 GB hard drive with at least 5 GB of available space

USB port or PCMCIA slot

CD or DVD-ROM drive

Windows-compatible monochrome or color printer

1024 x 768 or higher resolution video adapter and monitor

Ordering information	
EncorePro 2 and Data Cards	Part number
EncorePro 2 software	1054785
USB SmartCard reader/writer	1014352
PCMCIA SmartCard reader/writer	1014353
Encore SmartCards (pack of 10 cards)	1003304
SD card reader	1047300
SD card cover (attached to device)	1063858
SD cards (10 pack)	1063859
SD card mailer (10 pack)	1065146
Accessories	Part number
Encore reports manual	1055272 (Call customer service for other languages)
EncorePro 2 user guide	1060291
EncorePro 2 installation and administration guide	1054258
FOSQ 30 clinical reprint	1000141
FOSQ 30 questionnaire, English	1004550
Epworth Sleepiness Scale (pad of 50 sheets)	1058362

Respironics, Encore, and SmartCard are trademarks of Respironics, Inc. and its affiliates. All rights reserved. Microsoft, Windows, Windows Server, Windows Vista, and SQL Server are registered trademarks of Microsoft Corporation.



© 2009 Koninklijke Philips Electronics N.V.All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

 $\label{eq:CAUTION:US} \textbf{ Gederal law restricts these devices to sale by or on the order of a physician.}$

KW 8/18/09 MCI 4102719

Philips Healthcare is part of Royal Philips Electronics

Respironics Asia Pacific +65 6298 1088 Respironics Australia +61 (2) 9666 4444 Respironics Europe, Middle East, Africa +33 1 47 52 30 00 Respironics United Kingdom +44 800 1300 845 Philips Respironics 1010 Murry Ridge Lane Murrysville, PA 15668

Customer Service +1 724 387 4000 800 345 6443 (toll free, US only) www.philips.com/respironics