

A photograph of a male patient sitting at a table with his arm extended, while a female healthcare professional in blue scrubs stands behind him, adjusting a portable detector under his arm. In the background, there is a white mobile C-arm X-ray unit and a height scale on the wall.

PHILIPS

DuraDiagnost

Quality imaging, practical investment

As a total DXR solution provider, Philips has the right system for virtually every need. Our DuraDiagnost slips easily into the mix, delivering premium DR images in low acuity settings.

This floor-mounted, space saving system with a SkyPlate wireless portable detector is a sensible marriage of function and fiscal responsibility. High quality, high reliability facilitates confident diagnoses throughout the continuum of care.

Choose DuraDiagnost.

Key advantages

- Outstanding workflow efficiency
- High reliability for low cost of ownership
- Proven DR components deliver high-quality images
- Support from a global DR solutions provider

Reflecting the **need for quality** in low acuity settings


From every perspective, clinical, economic, imaging, and service, we've tailored DuraDiagnost to integrate seamlessly into its supporting role.

Clinical advantages

- Easy geometry positioning for outstanding workflow efficiency
- Intuitive to operate, the Eleva user interface helps establish a repeatable workflow, imperative for exam standardization
- Eleva uses patient data for quick and automatic pre-set exam parameters for exposure, image processing, export to PACS, and printing
- Accommodates a wide variety of chest, abdominal, and extremity exams
- UNIQUE image processing delivers consistently uniform clinical image quality

Economic value

- Space saving design for minimal room requirements
- Floor-mounted has lower installation costs compared to ceiling mounted systems
- Introduction of a digital workflow with rapid processing for enhanced efficiency
- SkyPlate wireless portable detector can be shared between compatible systems
- System uptime and reliability is enhanced with proven premium components, for a strong ROI
- Eleva user interface is easy to learn and simple to operate with only three clicks to get an image
- Because Eleva is common across Philips DXR systems, it offers the flexibility to use staff resources as necessary



Budgetary considerations may vary, but the need for quality never does. Designed to help you gain a competitive advantage and extend your clinical capacity, the DuraDiagnost offers an outstanding price to performance ratio.

When you acquire a Philips DuraDiagnost system, you can be certain it reflects the quality, reliability and history of a world-class brand.

- Images and information seamlessly integrate into most hospital networks
- DuraDiagnost reflects decades of company experience in digital X-ray development, and it is built for a low total cost of ownership

Imaging highlights

- Consistent, high quality digital images with the SkyPlate wireless detector mirror the imaging capabilities of Philips premium DR systems
- Lightweight SkyPlate wireless detector uses advanced Cesium Iodide (CsI) detector technology
- The Value Room uses the same quality tube and generator configuration as found on Philips premium DR systems



“With DuraDiagnost we are able to provide our patients with personalized medical care and at the same time remain within a reasonable budget.”

*Dr. Kathrin Scharfetter
Medical Specialist for Diagnostic
Radiology and Pediatric Radiology,
Radiology Center Eppendorfer Baum,
Hamburg, Germany*

Service solutions

- Take full advantage of Philips global services network
- Your RightFit service agreement is designed around your unique business priorities and capabilities
- With every RightFit service agreement get 24/7 access to our USA-based Customer Care Solutions Center staffed with a dedicated team of knowledgeable, responsive, and caring service associates
- Robust protection for wireless portable detector provides financial predictability
- Seamless OEM support for low total cost of ownership and predictable costs throughout the lifecycle of your equipment
- Optional Remote Services proactively monitors your systems for increased uptime

Quality system, affordable price

To facilitate clinical excellence and enhance ROI an X-ray system must be optimized for its specific role throughout the continuum of care. Low acuity environments for expanding healthcare networks, outpatient clinics, and overflow capacity require a quality system at an affordable price point. Philips DuraDiagnost fits the bill. Built with proven technology, it is designed for a low cost of ownership. Yet there is no image quality compromise. By using the very same digital imaging chain as our premium DR products, the DuraDiagnost Value Room delivers exceptional image quality, time after time.



Workflow features

- EasyAlign
- SmartOne button

Philips premium imaging chain

- Integrated generator
- Shares the same generator and tube as Philips premium DR

Floor mounted

- Durable design
- High uptime
- Low installation cost
- Low total cost of ownership

Introducing SkyPlate

- Large 14" x 17"
- SkyPlate sharing



New Eleva OS

- Based on WINDOWS 7
- UNIQUE
- Dose reporting (DAP)

© 2014 Koninklijke Philips N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.



Please visit www.philips.com/duradiagnost

Printed in The Netherlands.
4522 991 05661 * Aug 2014