



Philips AVENT  
Comfort Double electric  
breast pump

**Natural**

Includes 2x 4oz bottle



SCF334/02

## More comfort, more milk naturally

### Breast pump with massage cushion

When you are comfortable and relaxed, your milk flows more easily. That is why we created our most comfortable breast pump yet: Sit comfortably with no need to lean forward and let our soft massage cushion gently stimulate your milk flow.

#### **Sit more comfortably with no need to lean forward**

- More comfortable pumping position due to unique design

#### **More milk in less time**

- The perfect time saver for moms

#### **Choose your most comfortable setting**

- Features a gentle stimulation mode and 3 pumping settings

#### **Gently stimulates natural let down and milk flow**

- Soft massage cushion with warm feel

#### **Other benefits**

- Includes our Natural bottle and nipple for natural latch on
- Compact lightweight design, includes handy travel bag
- Intuitive assembly. Easy visual matching of parts
- Easy cleaning due to the small number of separate parts
- Easy operation at the touch of a button



asimpleswitch.com

**PHILIPS**  
**AVENT**

# Highlights

## More comfortable position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping; no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, naturally helps your milk to flow more easily.

## Simple settings to choose from



When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow most comfortable for you.

## Soft massage cushion



Our massage cushion has a new soft velvety texture that gives a warm feel to the skin for

comfortable, gentle stimulation of your milk flow. During pumping, it gently mimics your baby's suckling for natural let down and clinically proven effectiveness. The soft massage cushion is designed to comfortably fit the majority of moms. Should you need it, a larger size cushion is available separately (availability may vary per country).

## Easy cleaning

Cleaning is easy, thanks to the small number of separate parts. Your milk will never come into contact with the tubing and base unit. All parts are dishwasher proof, except the electrical parts.

## Easy operation

Double electric breast pumps are ideal for moms who pump milk regularly. Save time by comfortably pumping from both breasts at once.

## Perfect time saver for moms

Simultaneous pumping at both breasts is proven more efficient and may even boost your ability to produce breast milk.\*\*

## Intuitive assembly



Easy visual matching of parts for intuitive assembly.

## Compact design



The breast pump has a compact design, which makes it easy to hold and position on your breast. The small, lightweight base unit can easily be placed within comfortable reach for full control when pumping. For extra convenience in transport and storage, the tube simply wraps around the base unit. The breast pump comes with a handy travel bag.



[asimpleswitch.com](https://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Country of origin

- England

## Material

- Breast pump: Polypropylene, BPA free\*
- Bottle: Polypropylene, BPA free\*
- Nipple: Silicone, BPA free\*

## What is included

- Breast pump body: 2 pcs
- Base unit incl. tubing: 1 pcs
- Standard size cushion: 2 pcs
- Natural bottle 4oz: 2 pcs
- Extra soft Newborn Flow Nipple: 2 pcs
- Travel cover: 2 pcs
- Sealing disc for milk storage: 2 pcs
- Breast pad sample packs: 2 (2 Day pads & 2 Night pads) pcs
- Spare diaphragm: 2 pcs
- Travel bag: 1 pcs

## Breast pump

- Material: BPA free\*

## Design

- Breast pump design: Compact design
- Bottle design: Ergonomic shape, Wide neck

## Ease of use

- Breast pump use: Intuitive assembly, Easy cleaning, Fully compatible range

## Development stages

- Stages: Pregnancy, 0 - 6 months

## Functions

- No leaning forward: Sit in a comfortable position
- Settings: 3 Expression settings, 1 Stimulation mode
- More milk in less time: The perfect time saver
- Soft massage cushion: Gentle stimulation



Issue date 2014-12-17

Version: 5.0.1

© 2014 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* Following EU directive (2011/8/EU)

\*\* A randomised controlled trial to compare methods of milk  
expression after preterm delivery (Jones et al ADC 2001;85:F91).