





<b>PHILIPS</b>	<b>Product Recycling Passport</b> <b>LC0475 – WEEE Passport, FR3 Trainer</b>	Page 1 of 2
----------------	---	-------------


<b>Product name:</b>	<b>AED Trainer 3 (FR3 Defibrillator Trainer)</b>	
<b>Identification code(s)</b>	<b>861467, 989803159381</b>	
<b>Total weight (in Kg)</b>	<b>0.6 Kg</b>	
<b>Producer/ Manufacturer</b>	Name company:	Philips Medical Systems
	Address:	22100 Bothell-Everett Highway
	Zip code:	Bothell, WA 98021
	Country:	United States

<b>Recycle Info</b>	<b>Items:</b>	<b>Location</b>
<b>Special attention</b> 	Remove six (6) Phillips-head screws to open trainer.	<b>Figure 1</b>
	Contains circuit board.	<b>Figure 2</b>
<b>Fluids / Gases</b> 	<b>Items:</b>	<b>Location</b>
	N/A	
 <b>Batteries To be Removed</b>	<b>Type:</b>	<b>Location</b>
	Contains up to 4 AA alkaline batteries.	<b>Figure 1</b>
<b>Hazardous</b>  <b>To be Removed</b>	<b>Substances:</b>	<b>Location</b>
	N/A	


Remove AA Cells



**FIGURE 1**



Six Phillips  
head Screws



**FIGURE 2**

Remove Four Phillips head screws  
and remove Circuit Board